

Imaging Systems **FLUORESCENCE**

DOC-PRINT

> AS simple as:

1. switch on
2. press buttons
3. print or save

Applications

Nucleic acid detection

- Ethidium bromide
- SYBR™ Green

Protein detection

- Coomassie blue
- Silver Stain

Doc-Print is a stand-alone imaging system dedicated to basic gel documentation. Thanks to its built-in CompactFlash drive, the images can be saved on a CompactFlash or printed on a thermal printer for instant hard-copy images.

Doc-Print combines simplicity and economy for laboratories with budget in mind.

Key features

- No computer required
- Works with virtually any transilluminators
- Print and save in seconds
- Sensitive CCD camera
- Complimentary Photo-Capt software
- Extra large built-in LCD display for gel image preview



NEW
CompactFlash
drive



Real time



Integration
time



GLP



Electronic
Transilluminator



Free Software